

Bundle pluriStrainer® Mini 200 µm and Flow Cytometry Tubes



	Order No.	Size	Sterility	Packaging
☐	43-10200-46	25 pcs.	无菌型	25 pcs. (25 pcs. / Bag)
☐	43-10200-56	50 pcs.	无菌型	50 pcs. (25 pcs. / Bag)
☐	43-10200-66	100 pcs.	无菌型	100 pcs. (25 pcs. / Bag)
☐	43-10200-76	500 pcs.	无菌型	500 pcs. (25 pcs. / Bag)

描述

Bundle includes sterile pluriStrainers Mini with 100 µm mesh and sterile 5 ml polystyrene (PS) tubes (12x75mm). The combination of a flow cytometry tube and cell strainer offers a convenient way to prepare laboratory samples for flow cytometry analysis. A 5 µm polyethyleneterephthalat (PET) mesh is incorporated into the housing of the pluriStrainer Mini 5 µm. The polystyrene round-bottom tubes are closed with a cap. Tubes and the mini cell strainers will be delivered in separate bags with 25 pcs. each.

应用

- Obtaining a real single cell suspension after digestion of mammary tissue and organoids
- Preparation of a single cell suspension of blood cells from bone marrow, pancreas, thymus, lymph nodes and others
- Dissociation of cells from other primary tissue
- Preparation of real single cell suspension for flow cytometry (FACS™)

Additional Information

交付时间	1-2
网眼尺寸	200 µm
包装	25 pcs. (25 pcs. / Bag)
颜色	浅透明棕色
流量控制	否
Mesh Fixation	Injected in housing for strong hold
布料	PET(聚乙二醇酯)

外壳材料	PP（聚丙烯）
符合	1.5 ml reaction tube, 2.0 ml reaction tube, FACS™ tube, Cryo vials 5 ml / 10 ml, reaction tube with lid, 5 ml / 10 ml reaction tube with screw cap, 14- 15 ml conical centrifuge tube, 24 well plate, 48 well plate
处理	Gripping surface enables easy and aseptic handling.
不育症	无菌型
相当于	Cell Strainer Snap Cap™, Corning Cell Strainer, EASYstrainer™, Cell Sieve, Falcon™ Cell Strainer, Reversible™ Strainer, MACS® SmartStrainer, MACS® Pre-Separation Filters
稳定性	RLT® Buffer, TRIzol®, isopropanol, organic solvents
离心	是
组织解离	是
可消毒	是
运输条件	室温
储存条件	室温
法律信息	Buffer RLT® is a trademark of Qiagen, TRIzol® is a trademark of Molecular Research Center Inc., FACS™ is a trademark of BD Biosciences, Corning is a registered trademark of Corning Inc., EASYstrainer™ is a registered trademark of Greiner Bio One International GmbH, Cell Strainer Snap Cap™ and Falcon™ Cell Strainers is a registered trademark of Corning Inc., MACS® SmartStrainer is a registered trademark of Miltenyi Biotec GmbH

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.