

Monocytes Spin Medium


Order No.
Size

60-00095-10	100 ml
60-00095-11	250 ml
60-00095-12	500 ml

Description

Monocytes Spin Medium is a ready-to-use, sterile density medium for the enrichment of monocytes in a high yield from fresh peripheral blood or buffy coat through single step density centrifugation. Depending on the sample a purity of 60-80% is possible. Monocytes Spin is a gentle alternative to bead-mediated monocyte enrichments, especially for downstream applications like cell culture and differentiation.

Additional Information

Delivery Time (days)	1-2
Sample Material	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution
Storage Condition	Room temperature, protect from light
Regulatory Statement	For research use only. Not for use in diagnostic procedures.
pH Range	6.8 – 7.4

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.