



TDL - 40C

TDL - 50C

TDL - 60C

产品介绍

Production Introduction

该系列离心机壳采用高强度塑料，造型美观大方、高速运行噪音小、内腔温升小和分离效果好等优点。整机电脑控制，操作简便，配有电子门锁，安全可靠。使用于生物、医药学和农业科学等领域分析工作，在遗传基因，蛋白核算等实验研究的使用中效果更为显著。是医院、生物化学、实验室对血清血浆、放射免疫等作定性分析的理想仪器。

The machine`s frame is made of metal shell,so it`s weight is low,the noise is low when it ir running.The motro has a wide range of speed,and it also can run precisely.It can be operated simply and conveniently.It is safe and reliable because of being equipped with electronic lock.

It can be used for qualitative analysis ro serum,evivtion in the fields of biology,chemistry,medicine and so on.

主要特点

Main Features

- 无刷变频电机，力矩大，无粉尘、免维护。
 - 数字显示转速和时间。
 - 专利减震技术，运行平稳，噪音小。
 - 独特风道系统，内腔温升小。
 - 专利技术的自动门锁，安全可靠。
 - 专利技术的手动门锁，应急开门。
- Brushless frequency motor, torque, no dust, maintenance free.
 - Digital display speed and time.
 - Patent damping technology, stable running, low noise.
 - Unique air duct system, a small rise in the cavity temperature.
 - Automatic door lock - patent technology, safe and reliable.
 - The manual of patent technology open lock, emergency.



基本参数

型 号: Model	TDL-40C	TDL-50C	TDL-60C
定时范围: Timer Range	0min-99min	0min-99min	0min-99min
噪 音: Noise	< 70dB	< 70dB	< 70dB
电 源: Power	220V 50Hz	220V 50Hz	220V 50Hz
整机功率: Whole Power	400w	400w	400w
调速范围: Timing Range	100-4000rpm	100-5000rpm	100-6000rpm
外形尺寸: Dimensions	660×460×415mm (L×W×H)	595×400×360mm (L×W×H)	380×480×330mm (L×W×H)
重 量: Weight	70kg	60kg	60kg

1号角转子



2号角转子



3号角转子



配用转头

转子名称 Name	最高转速 Max speed	转子容量 Capacity	最大离心力 Max RCF	适用机型 Suitable model
1号角转子	4000rpm	15ml×24 50ml×8 100ml×4	3180×g	TDL-40C
2号角转子	5000rpm	15ml×8	3600×g	TDL-50C
3号角转子	6000rpm	15ml×12	5000×g	TDL-60C